

TURBIDITY CONVERSION TABLE

Purpose

To provide a lookup table to permit preparing bacterial inoculation samples for the Biolog Microbial Identification and Phenotype MicroArray™ products when the Biolog Turbidimeter is not available. The Biolog inoculation fluid is prepared in 20 mm OD test tubes. Cell suspensions are set up to specified densities using light scatter measurement. The suspended cells are added to the Biolog MicroPlates™ and the plates are incubated in a standard incubator or the OmniLog® Incubator/Reader.

GEN III

Conversion Table

Part Number	Turbidity Standard	Standard Value	Absorbance at 590 nm in a 1 cm cuvette
3411	GN-NENT	52 % T	0.128
3414	GP-COC & GP-ROD & GN-FAS	20 % T	0.315
3415	YT	47 % T	0.148
3419	GP-ROD SB	28 % T	0.249
3421	GN-ENT	61 % T	0.097
3426	FF	75 % T	0.056
3427	AN	65 % T	0.084
3431	PM: 85% T	85 % T	0.032
3432	PM: 40% T	40 % T	0.179

Microbial Identification	Inoculation Fluid Protocol	Cell Density Turbidity Range	Absorbance at 590 nm in a 1 cm cuvette
GEN III	A, B or C1	90% to 98%	0.020 to 0.004
GEN III	C2	62% to 68%	0.093 to 0.075

OmniLog is a registered trademark, MicroStation, MicroLog, MicroPlate, and RetroSpect are trademarks of Biolog, and all other trademarks are the property of their respective holders. Microsoft is a registered trademark of Microsoft, Inc.

TAS 013 Rev C

January 2012